Avicenum **ORTHO 1500**

(EN) INSTRUCTION MANUAL

Avicenum ORTHO 1500 Ankle Orthosis Type 01

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Avicenum ORTHO 1500 Sprunggelenkbandage Typ 01

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(EN) Dear customer,

In your hands, you are now holding an ankle orthosis for post-operative and post-traumatic therapy. You can maximize the benefits of the product by its proper use and maintenance.

Avicenum ORTHO 1500 Ankle Orthosis Type 01 is manufactured by the flat knitting technology. The materials used are non-irritating and latex-free. There are two anatomically shaped silicone pads integrated from the medial and lateral sides of the orthosis which improve proprioception (the ability to coordinate the body in space), thereby providing muscle support (ankle stabilization). The silicone pads have a massaging effect, helping to absorb swelling and blood clots. The mechanism behind the effects of the product is based on the compressive action in the ankle area which produces the desired therapeutic effect. In order to achieve the desired effects, you must select the correct size of the product. The product is most effective during physical strain. The product is also suitable for sports activities. You should remove the product during longer breaks.

Product name:

Avicenum ORTHO 1500 Ankle Orthosis Type 01

Product variants:

- Double-sided, silicone pads, sizes XS – XXL

User profile:

Avicenum ORTHO 1500 Ankle Orthosis Type 01 is a medical device intended for use both at home and in medical facilities.

Intended use:

A medical device designed to support the ankle during post-traumatic and post-operative therapy, in case of flabby muscles, inflammation, or injury to ankle ligaments and tendons and to prevent swelling. This product is also suitable for sports activities as a therapeutic or preventive aid.

Indications:

- Acute or chronic ankle instability
- Treatment after severe or mild distortions of the ankle (sprain)
- Ankle post-traumatic swelling
- Ankle swelling in arthritis or arthrosis
- Treatment after surgery or fracture of the ankle
- Painful syndromes of the ankle
- Prevention of ankle injury
- Ankle ligaments instability

If this is the first time you use the product, consult your doctor or have your medical condition regularly examined.

Contraindications:

- Skin diseases or injuries in the ankle area

- Acute weeping skin manifestations
- Hypersensitivity to materials used
- Lymphatic drainage disorders and unspecified swelling of the soft tissue in the ankle area
- Limb sensitivity disorders, peripheral neuropathy (e.g., in diabetes mellitus)

Your doctor should consider the use of the product in view of the known contraindications and possible risks compared to the product's therapeutic benefits.

Target group:

Patients with ankle instability, patients undergoing post-operative or post-traumatic therapy of the ankle, and patients with the above indications. The product is intended for use on healthy, intact skin. The product is intended for all age groups and both sexes. The product may be used in pregnancy without limitation.

Side effects:

Side effects may occur especially when using a wrong size of the product or if the product is used outside strenuous activities. If the product is too tight, there may be pressure sores, impaired blood circulation due to compression of blood vessels, or deterioration of sensitivity due to excessive compression of the nerves. Sensitive and hypersensitive individuals may experience skin irritation (redness, itching, or blistering) at the point of contact of the product with the skin.

How to choose the correct size:


In order to achieve the desired effects, you must select the correct size of the product. Determine the size by measuring the circumference at the points indicated in the figure.

Note:

Any swelling or instability of the ankle should be examined by your doctor. If this is the first time you experience swelling or instability of the ankle, consult your doctor before using the product. Do not use the product if you do not know the cause of the swelling or instability of the ankle. Always consult your doctor before combining the product with other medical devices (e.g., compression stockings).

- Use a tape measure.
- Take the measurements standing up, preferably in the morning right after waking up.
- Get assistance from another person.
- Measured values may change during the therapy – always verify the correct size before using a new product.

Measured above the ankle		
XS	S	M
19 - 21	21 - 23	23 - 25
L	XL	XXL
25 - 27	27 - 29	29 - 31

If you are on the border between two sizes, choose the bigger one. If you have difficulties choosing the size of the product, visit a specialized medical supplies store or directly the manufacturer.

If you have difficulties putting on the product, visit a specialized medical supplies store and request a demonstration of the correct method of applying it. You can also use one of the aids for putting on the product.

With a correctly selected size:

- You should feel noticeable but not unpleasant pressure in the ankle area
- The product does not slide or roll down when walking
- The product does not choke at any point
-

Size and manufacturer information can be found on the tag sewn at the top edge of the product

Method of use:

Put on the product by gradually pulling it over the ankle so that the silicon supports fit to both the lateral and medial malleoli. Never put on the orthosis by merely pulling the top edge – always put the orthosis on gradually. It is best to put on the orthosis in the morning (while sitting) before any physical strain. When putting on, the product may not come into contact with sharp objects (e.g., jewellery) or long nails. Intensive care for your feet (no rough skin on the heels) will decrease the risk of damaging the product and increase its life. Stretch the orthosis evenly without creases for effective distribution of pressure.

If you use body cosmetics, wait about 15 minutes after application before you put on the orthosis.

In the case of visible mechanical damage to the orthosis, stop wearing it and preferably replace it with a new product.

Care and maintenance:

Avicenum ORTHO 1500 Ankle Orthosis Type 01 is a medical device intended for repeated use by a single person. The product will maintain its properties and safety for 12 months from the first use under these conditions:

- The product can be washed in a washing machine (maximum washing temperature 30 °C, gentle cycle) in a soapy solution or using a special preparation for washing compressing stockings and orthosis using a wash bag
- Do not use fabric softeners
- Do not use centrifugation; squeeze the wet product between two towels
- Dry the product in a horizontal position away from a direct heat source (radiator, sun)
- Do not iron the product.
- The maximum number of washing cycles is 120
- Store in a dry and dark place, preferably in the original packaging
- Do not interfere mechanically in the product, do not repair places damaged by use, do not cut the edges of the product
- The product may not come into contact with organic solvents or bleaching agents

The medical device may be used safely for 24 months from the date of manufacture. This date is contained in the batch number of the device expressed as YY/MM-1234567 where YY is the year and MM the month of manufacture.

Care symbols:

Do not use fabric softeners

Prohibited use:

Stop using the product immediately if:

- If it shows signs of mechanical damage (replace the product)
- Your health has deteriorated which could be associated with the use of the product – consult your doctor at once

Composition:

60 % PAD Nylon, 21 % elastane LYCRA®, 19 % PES Polyester, cotton terry, silicone

This medical device does not contain an active substance, including human blood or plasma derivatives, tissues or cells of human origin or derivatives thereof, tissues or cells of animal origin or derivatives thereof.

Disposal:

The product can be disposed of with normal household waste or in a textile waste container. Although the product is intended for use on healthy and intact skin, it may be contaminated with body fluids during use. In that case, dispose of the product as an infectious material.

The material used is not toxic unless burned. Do not burn the used orthosis.

Any serious adverse event that occurs in connection with the device should be reported to the manufacturer and the competent authority of the Member State in which the user/patient is resident.

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Manufacturer:

ARIES, a.s., Studenec 309, 512 33 Studenec
Last reviewed: 01.01.2020

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