Avicenum **ORTHO 2000**

(EN) INSTRUCTION MANUAL
AVICENUM ORTHO 2000 Wrist orthosis type 01

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(EN) Dear customer,

you are holding a medical wrist orthosis to support post-operative and post-traumatic therapy. You can maximise the benefits of this product by its proper use and maintenance.

Avicenum ORTHO 2000 wrist orthosis Type 01 is made of 3D fabric. The materials used for its production are non-irritating and latex-free. In the palm part of the orthosis, there is an anatomically shaped, removable splint and two other smaller splints on the sides of the wrist, which in combination with the hook-and-loop fasteners improve support and fixation of the wrist joint. The mechanism of action of the wrist orthosis is based on the compression action and fixation in the wrist area, which produces the necessary therapeutic effect. A proper fit is necessary to achieve the desired effect. For longer rest breaks, you should remove the orthosis unless otherwise directed by your attending doctor.

Product name:

AVICENUM ORTHO 2000 wrist orthosis Type 01

Product variants:

- left, size S–XL
- right, size S–XL

User profile:

Avicenum ORTHO 2000 wrist orthosis Type 01 is a medical device intended for use both at home and in medical care facilities.

Intended use of the medical device:

A medical device designed to support the wrist during post-traumatic and post-operative therapy, prevent swelling and serve as support for wrist ligament and tendon laxity, inflammation or injury.

Indications:

- acute or chronic instability of the wrist
- follow-up treatment after severe or less severe wrist distortions (sprains)
- post-traumatic swelling of the wrist
- swelling of the wrist in arthritis or arthrosis
- enthesopathy in the radioulnar articulation
- follow-up treatment after surgeries and fractures of the wrist
- painful wrist syndromes
- tendonitis in the wrist area
- carpal tunnel syndrome
- prevention of wrist injuries
- instability of the ligamentous apparatus of the wrist

If you are using the wrist orthosis for the first time, it is advisable to consult your attending doctor or have your health checked regularly by your doctor.

Contraindications:

- skin diseases or injuries in the wrist area
- acute weeping skin manifestations
- hypersensitivity to materials used
- lymphatic drainage disorders and non-specific soft tissue swelling in the wrist area
- limb numbness disorders, peripheral neuropathy

In view of known contraindications your attending doctor should evaluate possible risks and weigh them against the therapeutic benefits of use of medical device when considering application of the wrist orthosis.

Target group:

Patients with wrist instability, post-operative and post-traumatic therapy of the wrist and patients with the above mentioned indications. The wrist orthosis is intended for use on healthy, intact skin. The orthosis is designed for all age groups and both genders. There are no restrictions for use during pregnancy.

Side effects:

Side effects may occur, especially if the wrist orthosis does not fit properly or if it is worn when at rest. There may be bruising, impaired blood flow due to pressure on blood vessels, or impaired sensitivity due to a too tightly fitting wrist orthosis. Sensitive and hypersensitive individuals may experience skin irritation (redness, itching, or blistering) at the place of contact of the wrist orthosis with the skin.

How to choose the correct size of the wrist orthosis:


A proper fit is necessary to achieve the desired effect. Determine the size by measuring the circumference of the limb at the points indicated in the figure.

Note:

Any swelling or instability of the wrist should be assessed by a doctor. If this is the first time you have noticed swelling or instability of the wrist, consult your doctor before using the wrist orthosis. Do not start self-treatment with the wrist orthosis unless you know the cause of the swelling or instability of the wrist. Always consult your attending doctor before using this product in combination with other medical devices (such as compression arm sleeve).

- use a tape measure
- ask someone to help you
- measured values may change during the therapy – always verify the correct size before using a new wrist orthosis

Measured wrist circumference point "c"	
S	M
12 - 16	16 - 20
L	XL
20 - 24	24 - 28

If you are on the border between two sizes, choose the bigger one. If you are unsure of your size of your wrist orthosis, visit a specialised medical supplies store or consult the manufacturer.

If you have difficulties putting on the wrist orthosis, visit a specialised medical supplies store and request a demonstration of the correct method of applying it.

With the correctly selected size:

- you feel a noticeable but not uncomfortable pressure in the wrist area
- the wrist orthosis does not slip or come down
- the wrist orthosis does not pinch at any point

Size and manufacturer information can be found on the label sewn to the top edge of the wrist orthosis.

Application of the wrist orthosis:

Push your thumb through the hole in the orthosis. The orthosis is placed on the wrist so that the support splints sit on the wrist. Tighten the orthosis using the hook-and-loop fasteners so that the pressure on the wrist is not uncomfortable. When putting on the orthosis, it should not come into contact with sharp objects (e.g. jewellery) or long nails.

The orthosis must be stretched evenly, without folds, so as to achieve the efficient pressure distribution. If you use body cosmetics, wait about 15 minutes after application before putting on the orthosis.

Should the orthosis become visibly mechanically damaged, stop wearing it and replace it with a new one.

Care and maintenance:



Avicenum ORTHO 2000 wrist orthosis Type 01 is a medical device intended for repeated use by one person. The product will maintain its properties and safety for 12 months from the first use under these conditions:

- the wrist orthosis can be machine washed (after removing the splints in the palm part of the orthosis) (max. temperature 30 °C, gentle cycle) with a soapy solution or a special detergent for compression stockings and bandages, in a protective cover (washing bag)
- always leave the hook-and-loop fasteners fastened when washing
- do not use fabric softener
- do not use centrifugation, squeeze the wet wrist orthosis between two towels

- lay flat to dry, away from a direct heat source (radiator, direct sunlight)
- do not iron the wrist orthosis
- maximum number of washing cycles is 60
- store in a dark, dry place in the original packaging
- do not alter the wrist orthosis in any way, do not repair areas damaged by use, do not trim the edges
- keep the wrist orthosis away from organic solvents and bleaching agents

This medical device may be used safely for 36 months from the date of manufacture. This date is part of the product batch number expressed in the format YY/MM-1234567, where YY is the year and MM the month of manufacture.

Care symbols:

When not to use:

Stop using the medical device immediately if any of the following conditions apply:

- it shows sign of mechanical damage (replace the orthosis with a new one)
- you notice a deterioration in your health that may be associated with the use of the wrist orthosis – consult your attending doctor in such case

Composition:

PAD Nylon, elastane, PES, metal, plastic

This medical device does not contain any active substance, including human blood or plasma derivatives, tissues or cells of human origin or derivatives thereof, or tissues or cells of animal origin or derivatives thereof.

Disposal:

The wrist orthosis can be disposed of with normal household waste or recycled as a textile. Although the product is intended for use on healthy, intact skin, it may become contaminated with body fluids during use. In this case, the wrist orthosis should be disposed of as infectious material.

The material used is non-toxic unless ignited. Never burn a used orthosis.

Any serious adverse event occurring in connection with the product in question should be reported to the manufacturer and to the competent authority of the Member State where the user and/or patient is located:

Czech Republic	State Institute for Drug Control, Šrobárova 48, 100 41 Prague 10
Slovakia	State Institute for Drug Control, Kvetná 1024/11, 821 08 Bratislava
Poland	Office of Registration of Medical Products, Medical Products and Biological Products, Al. Jerozolimskie 181C, 02-222 Warsaw
Germany	Federal Institute for Drugs and Medical Devices, Kurt Georg Kiesinger Allee 3, 53175 Bonn
Russia	Roszdraznadzor, Address: 4, bld. 1, Slavyanskaya Square, Moscow, 109074
Ukraine	Ministry of Health Ukraine, 7 Hrushevskoho Street, Kyiv, 01601
Australia	Therapeutic Goods Administration, 136 Narrabundah Lane, Symonston ACT 2609

Manufacturer:

ARIES, a.s., Studenec 309, 512 33 Studenec
 Last reviewed: 17.05.2021

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