Avicenum **140**

(GB) Instructions for Use

Avicenum PHLEBO 140 FINE knee-high socks
Avicenum PHLEBO 140 STRONG knee-high socks
Avicenum PHLEBO 140 FINE Thigh-High Stockings
Avicenum PHLEBO 140 STRONG Thigh-High Stockings
Avicenum PHLEBO 140 FINE Tights

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Avicenum PHLEBO 140 FEINE Wadenstrümpfe
Avicenum PHLEBO 140 DICKE Wadenstrümpfe
Avicenum PHLEBO 140 FINE Schenkelstrümpfe
Avicenum PHLEBO 140 DICKE Schenkelstrümpfe
Avicenum PHLEBO 140 FEINE Strumpfhose

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(GB)

INSTRUCTIONS FOR USE

The medical compression stockings that you hold in your hands are a useful device in the treatment of circulatory and lymphatic disorders. You can maximize the benefits of this product by its proper use and maintenance.

Avicenum PHLEBO 140 stockings are a compression class 1 medical device (I. KT, CCL 1). They were knitted and do not contain latex. This compression medical device works by applying precisely defined, made by circular knitting technology. The materials used in our compression medical devices are kind to the graduated pressure to the leg, with a gradual reduction in force from the ankle towards the heart. The pressure is greatest in the ankle area at point b (see size chart), reaching values of 18-21 mmHg. A proper fit is necessary to achieve the desired effect. Use the compression stocking when engaged in physical activity; compression products should be removed for longer rest periods (unless otherwise indicated by your doctor). **Avicenum PHLEBO 140** is a medical device intended for use both at home and in medical care facilities.

Product name:

Avicenum PHLEBO 140 FINE Knee-high socks

Product variants:

- Closed toe, I.KT, size S normal – L long

Avicenum PHLEBO 140 STRONG Knee-high socks

Product variants:

- Closed toe, I.KT, size S normal – L long

Avicenum PHLEBO 140 FINE Thigh-high stockings

Product variants:

- Lace, closed toe, I.KT, size S normal – L long

Avicenum PHLEBO 140 STRONG Thigh-high stockings

Product variants:

- Hem, closed toe, I.KT, size S normal – L long

Avicenum PHLEBO 140 FINE Tights

Product variants:

- Closed toe, I.KT, size S normal – L long

User profile:

Avicenum PHLEBO 140 is a medical device intended for use both at home and in medical care facilities.

Intended use of the medical device:

This medical device is intended for compression therapy of circulatory disorders and various types of oedema in the leg.

Indications:

- Chronic venous insufficiency – its subjective and objective manifestations, prevention
- Primary and secondary varicose veins in the legs
- Superficial thrombophlebitis
- Oedemas associated with post-thrombotic syndrome after overcoming phlebothrombosis
- Oedemas and varicose veins in pregnancy
- Conditions following the occurrence of phlebothrombosis and thrombophlebitis in the legs
- Prevention of phlebothrombosis in surgical treatment and subsequently also thromboembolic conditions (especially surgery, orthopaedics, gynaecology)
- Oedemas of the lower limbs (postoperative and post-traumatic, lymphatic, obesity-related oedemas, idiopathic oedemas of the lower limbs)
- Compression support following sclerotherapy of varicose veins of the legs
- Compression support following surgical extirpation of varicose veins (classic surgery, radiofrequency ablation or laser ablation of varicose veins)
- Travel-related thrombosis (economy class syndrome)

If you are using a compression medical device for the first time, consult your doctor to determine the proper compression class and have regular check-ups during the treatment.

Contraindications:

- Acute weeping skin manifestations
- Acute limb ischaemia
- Critical limb ischaemia
- Decompensated cardiac insufficiency
- Phlegmasia cerulea dolens, phlegmasia alba dolens
- Hypersensitivity to materials used
- Limb numbness disorders (e.g. peripheral neuropathy in diabetes mellitus)

Your doctor should consider the use of the compression stockings in view of the known contraindications and possible risks compared to the product's therapeutic benefits.

Target group:

Patients with primary varicose veins, patients after venous thrombosis and patients with the above-mentioned indications. Compression medical stockings are intended for use on healthy, intact skin. This compression medical device is suitable for all age groups and both sexes. There are no restrictions for use during pregnancy.

Side effects:

Side effects may occur, especially if the compression medical device does not fit properly or if it is worn when at rest. There may be bruising, impaired blood flow due to pressure on blood vessels, or impaired sensitivity due to excessive nerve compression. Sensitive and hypersensitive individuals may experience skin irritation

(redness, itching, or blistering) in places where the compression stocking comes into contact with the skin.

Fitting your compression stockings:

A proper fit is necessary to achieve the desired effect. Use the size chart to determine the correct size by measuring limb circumference at the points indicated.

Note:

Any oedema of the lower limbs should be examined by your doctor. If this is the first time you experience oedema, consult your doctor before using the product. Do not use the product if you have not determined the cause of the oedema. Always consult your doctor before using this product in combination with other medical devices (such as bandages or orthoses).

- Use a tape measure
- Take measurements while standing, preferably first thing in the morning – legs may swell over the course of the day
- Ask someone to help you
- Measurements may change over the course of treatment – always check your size before using a new pack of compression stockings


SIZE TABLE:

SIZE	b	d	f	g	h	t
S	20 – 24	34 – 37	46 – 50	50 – 55	110	80
M	24 – 27	37 – 40	50 – 54	55 – 60	120	90
L	27 – 30	40 – 44	54 – 58	60 – 66	130	100

LENGTH			
	short	normal	long
AD	–	34 – 38	38 – 43
AG	56 – 62	62 – 71	71 – 82

- a foot circumference
 b circumference above the ankle
 b1 circumference 10 cm above the ankle
 c calf circumference
 d circumference below the knee
 e circumference above the knee
 f thigh circumference (in the middle of the thigh)
 g thigh circumference (5 cm under the crotch)
 h hip circumference
 t waist circumference

- A – D knee-high socks
 A – G thigh-high stockings
 A – T thighs

If the measured values do not correspond to any of the sizes in the table (e.g. point b falls under size S, but point d is size L), visit a specialised medical supplies store or consult the manufacturer. You may need a custom-made compression stocking.

If the measured values differ significantly for the right and left limb, visit a specialised medical supplies store or consult the manufacturer. You may need a custom-made compression stocking.

If you are on the border between two sizes, choose the bigger one. If you are unsure of your size, visit a specialised medical supplies store or consult the manufacturer.

If you have difficulty putting on your compression medical stockings, visit a specialised medical supplies store and request a demonstration of the correct method of putting it on. You can also use one of the aids for putting on the product.

If the hem or lace does not stay in place on the leg (for thigh-high stockings), wipe the silicone strip or studs with an alcohol solution, or shave the hair on the limb at the point of contact of the silicone with the skin.

With the correctly selected size:

- You should feel noticeable but not unpleasant pressure over the entire length of the compression medical device
- The compression medical device should not slide down when walking and should not bunch up under the knee or at the instep
- The hem or lace edge of the compression stockings should not be too tight (for thigh-high stockings)

Information on compression class, size, shelf life and manufacturer can be found on the label on the medical device packaging.

Putting on your compression medical stockings:

Pull the compression stocking onto the leg gradually; do not tug on the upper edge alone. Compression medical stockings should be put on in the morning before any physical activity and before any oedema occurs. When pulling on the compression medical stockings, keep them away from sharp objects (e.g. jewellery) or long nails. Increased care of the soles of your feet (rough skin on the heels) will reduce the risk of damaging the product. The compression stockings must be stretched evenly, without folding, to ensure effective pressure distribution. If you use any body cosmetics, allow 15 minutes after application before putting on the medical device. In the event of visible mechanical damage to the product, stop wearing it and replace it with a new one.

Care and maintenance:

Avicenum PHLEBO 140 compression medical stockings are intended for repeated use by a single individual. This product may be used safely and effectively for 6 months from the first use if the following guidelines are respected:

- This compression medical device can be machine washed (max. temperature 30°C, gentle cycle) with a soapy solution or a special detergent for compression medical devices, in a protective cover (washing bag)
- Do not use fabric softener
- Do not spin-dry, press the wet compression medical device dry between two towels
- Lay flat to dry, away from a direct heat source (radiator, direct sunlight)
- Do not iron compression medical devices
- Maximum number of washing cycles is 120
- Store in a dark, dry place in the original packaging
- Do not alter the compression medical device in any way, do not repair areas damaged by use, do not trim the lace or hem
- Restore the adhesive quality of the silicone strips on the lace or the silicone nubs on the hem by wiping with an alcohol solution
- Keep the compression medical device away from organic solvents and bleaching agents

This medical device may be used safely for 60 months from the date of manufacture. This date is part of the product batch number expressed in the format YY/MM-1234567, where YY is the year and MM the month of manufacture.

Care symbols:

Do not use fabric softeners!

When not to use:

Stop using the product immediately if any of the following conditions apply:

- There are signs of mechanical damage – replace with new compression medical stockings
- You notice a deterioration in your health that may be associated with the use of the compression medical device – consult your doctor in such case

Composition:

Knee-highs stockings: 69% PAD nylon, 31% elastane LYCRA®

Thigh-high stockings: 69% PAD nylon, 31% elastane LYCRA®, silicone

Tights: 69% PAD nylon, 31% elastane LYCRA®

This medical device does not contain any active substance, including human blood or plasma derivatives, tissues or cells of human origin or derivatives thereof, or tissues or cells of animal origin or derivatives thereof.

Disposal:

Compression medical stockings can be disposed of with normal household waste or in a textile container. Although the product is intended for use on healthy, intact skin, it may become contaminated with body fluids during use. In such a case, dispose of the compression medical device as infectious material.

The material used is non-toxic unless ignited. Never incinerate used compression medical devices.

Any serious adverse event occurring in connection with the product in question should be reported to the manufacturer and to the competent authority of the Member State where the user and/or patient is located:

Czech Republic	Státní ústav pro kontrolu léčiv, Šrobárova 48, 100 41 Prague 10
Slovak Republic	Štátny ústav pre kontrolu liečiv, Kvetná 1024/11, 821 08 Bratislava
Poland	Urząd Rejestracji Produktów Leczniczych, Wyrobów Medycznych i Produktów Biobójczych, Al. Jerozolimskie 181C, 02-222 Warsaw
Germany	Das Bundesinstitut für Arzneimittel und Medizinprodukte, Kurt Georg Kiesinger Allee 3, 53175 Bonn
Russia	Roszdraznadzor, Address: 4, bld. 1, Slavyanskaya Square, Moscow, 109074
Ukraine	Ministry of Health Ukraine, 7 Hrushevskoho Street, Kiev, 01601
Australia	Therapeutic Goods Administration, 136 Narrabundah Lane, Symonston ACT 2609

List of abbreviations:

KT, CCL – compression class

mmHg – millimeters of mercury column (unit of blood pressure measurement on the leg)

Manufacturer:

ARIES, a.s., Studenec 309, 512 33 Studenec
Last reviewed: 31.10.2022

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Ukrajina	Ministry of Health Ukraine, 7 Hrushevskoho Street, Kyiv, 01601
Austrália	The Therapeutic Goods Administration, 136 Narrabundah Lane, Symonston ACT 2609

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We proudly certify that

Aries A.S.

has certified all medical devices that the respective articles have successfully passed all tests according to the stringent "Aries" Aries Standard and that it is certified as a regular quality control program.

The company has been approved as "Supplier" license partner and has the right to distinguish its respective articles with the "Aries" quality seal.

Authorized licensee of the Aries Standard certificate

OEKO-TEX®
STANDARD 100
by OEKO-TEX®

OEKO-TEX®
CERTIFIED BY TESTEX

CERTIFICATE

The company

Aries, a.s.
15200 Studenec uvalna, 02038 NEPALSKO

is producing/certifying according to STANDARD 100 by OEKO-TEX® for use in CATEGORY: 001 by OEKO-TEX® (Textiles, leather or fur apparel) 0000-00002

OEKO-TEX®
STANDARD 100
by OEKO-TEX®

for the following articles:

Medical equipment and supporting goods (linens, anti-static), combination sport products, sport and leisure clothing, outdoor gear in the form of, for and with linens, which prove freedom of allergenic and/or genetic, according to STANDARD 100-000, and/or that include already certified according to standard 001 by OEKO-TEX®

Products in the respective main categories (0000001-001 to 0000-00002), however it is possible that they are also produced by other OEKO-TEX® partner manufacturing companies (OEKO-TEX® 0000001-0000000) - OEKO partner's responsibility regarding the product will then remain in place.

This certified product falls respectively into Category 001 of STANDARD 100. Use of non-certified, other articles may also utilize the Aries or Aries partner's responsibility system of use of Aries Standard (Aries) with the exception of the above-mentioned items.

Issuance of this certificate, which has been a condition of certification according to STANDARD 100-000, is subject to an obligation from the licensee (Aries) to OEKO-TEX® and only in cooperation with projects for sustainable development which relate to the industry to which it refers.

This certificate 000001 is valid until 31.03.2018

Since 15.03.2017

"OEKO-TEX" registered for OEKO-TEX, Vienna, Austria.

Aries A.S. Studenec uvalna 309, 15200 Studenec, NEPALSKO

Aries
MEDICAL PRODUCTS

ARIES, a.s., 512 33 Studenec 309
MADE IN CZECH REPUBLIC
ISO 9001:2015

www.avicenum.eu