Avicenum **PHLEBO 250**

(GB) Instructions for Use

Avicenum PHLEBO 250 Calf-Length Stockings
Avicenum PHLEBO 250 Thigh-High Stockings
Avicenum PHLEBO 250 Tights

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Avicenum PHLEBO 250 Wadenstrümpfe
Avicenum PHLEBO 250 Schenkelstrümpfe
Avicenum PHLEBO 250 Strumpfhose

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(GB)

INSTRUCTION MANUAL

Dear customer, the medical compression stockings that you hold in your hands are a useful aid in the treatment of circulatory and lymphatic disorders. You can maximize the benefits of this product by its proper use and maintenance.

Avicenum PHLEBO 250 stockings are a compression class 1 medical aid (I.KT, CCL1). They were made by circular knitting technology. The materials used in our compression medical aids are kind to the skin and do not contain latex. This compression medical aid works by applying precisely defined, graduated pressure to the leg, with a gradual reduction in force from the ankle towards the heart. The pressure is greatest in the ankle area at point b (see size chart), reaching values of 18 - 21 mmHg. A proper fit is necessary to achieve the desired effect. Use the compression stocking when engaged in physical activity; compression products should be removed for longer rest periods (unless otherwise indicated by your doctor). **Avicenum PHLEBO 250** is a medical aid intended for use both at home and in medical care facilities.

Product name:

Avicenum PHLEBO 250 Calf-Length Stockings

Product variants:

- Open toe, I.KT, size S normal – XL long
- Closed toe, I.KT, size S normal – XL long

Avicenum PHLEBO 250 Thigh-High Stockings

Product variants:

- Lace, open toe, I.KT, size S normal – XL long
- Lace, closed toe, I.KT, size S normal – XL long
- Plain hem, open toe, I.KT, size S normal – XL long
- Plain hem, closed toe, I.KT, size S normal – XL long

Avicenum PHLEBO 250 Tights

Product variants:

- Closed toe, I.KT, size S normal – XL long

Product accessories:

Each package of **Avicenum PHLEBO 250** medical stockings, including all variants, contains a textile dressing aid for Avicenum compression stockings. The dressing aid is not a medical device and serves exclusively to assist you in putting on your Avicenum compression stockings. The use of the dressing aid is described below, along with instructions for how best to put on compression medical devices.

User profile:

Avicenum PHLEBO 250 is a medical aid intended for use both at home and in medical care facilities.

Intended use:

This medical aid is intended for compression therapy of circulatory disorders and various types of oedema in the leg.

Indications:

- chronic venous insufficiency – its subjective and objective manifestations, prevention
- primary and secondary varicose veins in the legs
- superficial thrombophlebitis
- oedemas associated with post-thrombotic syndrome after overcoming phlebothrombosis
- oedemas and varicose veins in pregnancy
- conditions following the occurrence of phlebothrombosis and thrombophlebitis in the legs
- prevention of phlebothrombosis in surgical treatment and subsequently also thromboembolic conditions (especially surgery, orthopaedics, gynaecology)
- oedemas of the lower limbs (postoperative and post-traumatic, lymphatic, obesity-related oedemas, idiopathic oedemas of the lower limbs)
- compression support following sclerotherapy of varicose veins of the legs
- compression support following surgical extirpation of varicose veins (classic surgery, radiofrequency ablation or laser ablation of varicose veins)
- travel-related thrombosis (economy class syndrome)

If you are using a compression medical aid for the first time, consult your doctor to determine the proper compression class and have regular check-ups during the treatment.

Contraindications:

- acute weeping skin manifestations
- acute limb ischaemia
- critical limb ischaemia
- decompensated cardiac insufficiency
- phlegmasia cerulea dolens, phlegmasia alba dolens
- hypersensitivity to materials used
- limb numbness disorders (e.g. peripheral neuropathy in diabetes mellitus)

Your doctor should consider the use of the compression stockings in view of the known contraindications and possible risks compared to the product's therapeutic benefits.

Target group:

Patients with primary varicose veins, patients after venous thrombosis and patients with the above-mentioned indications. Compression medical stockings are intended for use on healthy, intact skin. This compression medical aid is suitable for all age groups and both sexes. There are no restrictions for use during pregnancy.

Side effects:

Side effects may occur, especially if the compression medical aid does not fit properly or if it is worn when at rest. There may be bruising, impaired blood flow due to pressure on blood vessels, or impaired sensitivity due to excessive nerve compression. Sensitive and hypersensitive individuals may experience skin irritation (redness, itching, or blistering) in places where the compression stocking comes into contact with the skin.

Fitting your compression stockings:

A proper fit is necessary to achieve the desired effect. Use the size chart to determine the correct size by measuring limb circumference at the points indicated.

Note:

Any oedema of the lower limbs should be examined by your doctor. If this is the first time you experience oedema, consult your doctor before using the product. Do not use the product if you have not determined the cause of the oedema. Always consult your doctor before using this product in combination with other medical aids (such as bandages or orthoses).

- use a tape measure
- take measurements while standing, preferably first thing in the morning – legs may swell over the course of the day
- ask someone to help you
- measurements may change over the course of treatment – always check your size before using a new pack of compression stockings


SIZE TABLE:

SIZE	a	b	b1	c	d	e	f	g	h	t
S	19 – 25	20 – 23	25 – 30	32 – 38	31 – 37	33 – 39	41 – 49	47 – 57	110	80
M	22 – 28	23 – 26	28 – 33	35 – 41	34 – 40	36 – 42	45 – 53	52 – 62	120	90
L	25 – 31	26 – 29	31 – 36	38 – 44	37 – 43	39 – 45	49 – 57	57 – 67	120	90
XL	28 – 34	29 – 32	34 – 39	41 – 47	40 – 46	42 – 48	53 – 61	62 – 72	130	100

LENGTH		
	normal	long
AD	34 – 38	38 – 43
AG	62 – 71	71 – 82

- a foot circumference
- b circumference above the ankle
- b1 circumference 10 cm above the ankle
- c calf circumference
- d circumference below the knee
- e circumference above the knee
- f thigh circumference (in the middle of the thigh)
- g thigh circumference (5 cm under the crotch)
- h hip circumference
- t waist circumference

- AD calf-length stockings
- AG thigh-high stockings

If the measured values do not correspond to any of the sizes in the table (e.g. point b falls under size S, but point d is size L), visit a specialised medical supplies store or consult the manufacturer. You may need a custom-made compression stocking.

If the measured values differ significantly for the right and left limb, visit a specialised medical supplies store or consult the manufacturer. You may need a custom-made compression stocking.

If you are on the border between two sizes, choose the bigger one. If you are unsure of your size, visit a specialised medical supplies store or consult the manufacturer.

If you have difficulty putting on your compression medical stockings, visit a specialised medical supplies store and request a demonstration of the correct method of putting it on. You can also use one of the aids for putting on the product.

If the hem or lace does not stay in place on the leg (for thigh-high stockings), wipe the silicone strip or nubs with an alcohol solution, or shave the hair on the limb at the point of contact of the silicone with the skin.

With a correctly selected size:

- you should feel noticeable but not unpleasant pressure over the entire length of the compression medical aid
- the compression medical aid should not slide down when walking and should not bunch up under the knee or at the instep
- the plain or lace edge of the compression stockings should not be too tight (for thigh-high stockings)

Information about compression class, size, and manufacturer can be found on the tag sewn at the upper edge of the stocking. Should the tag be removed in the course of use, sizing information is also knit in on the stocking sole.

Putting on your compression medical stockings:

Pull the compression stocking onto the leg gradually; do not tug on the upper edge alone. Compression medical stockings should be put on in the morning before any physical activity and before any oedema occurs. For easier application, use the enclosed Avicenum compression stocking textile dressing aid. The use of the textile dressing aid is illustrated on the enclosed card. When pulling on the compression medical stockings, keep them away from sharp objects (e.g. jewellery) or long nails. Increased care of the soles of your feet (rough skin on the heels) will reduce the risk of damaging the product. The compression stockings must be stretched evenly, without folding, to ensure effective pressure distribution. If you use any body cosmetics, allow 15 minutes after application before putting on the medical aid. In the event of visible mechanical damage to the product, stop wearing it and replace it with a new one.

Care and maintenance:

Avicenum PHLEBO 250 compression medical stockings are intended for repeated use by a single individual. This product may be used safely and effectively for 6 months from the first use if the following guidelines are respected:

- this compression medical aid can be machine washed (max. temperature 30 °C, gentle cycle) with a soapy solution or a special detergent for compression medical aids, in a protective cover (washing bag)
- do not use fabric softener
- do not spin-dry, press the wet compression medical aid dry between two towels
- lay flat to dry, away from a direct heat source (radiator, direct sunlight)
- do not iron compression medical aids
- maximum number of washing cycles is 120
- store in a dark, dry place in the original packaging
- do not alter the compression medical aid in any way, do not repair areas damaged by use, do not trim the lace or edge
- restore the adhesive quality of the silicone strips on the lace or the silicone nubs on the hem by wiping with an alcohol solution
- keep the compression medical aid away from organic solvents and bleaching agents

This medical aid may be used safely for 54 months from the date of manufacture. This date is part of the product batch number expressed in the format YY/MM-1234567, where YY is the year and MM the month of manufacture.

Care symbols:

Do not use fabric softeners!

When not to use:

Stop using the product immediately if any of the following conditions apply:

- there are signs of mechanical damage – replace with new compression medical stockings
- you notice a deterioration in your health that may be associated with the use of the compression medical aid – consult your doctor in such case

Composition:

Calf-length stockings: 65 % PAD nylon, 35 % elastane LYCRA®

Thigh-high stockings: 65 % PAD nylon, 35 % elastane LYCRA®, silicone

Stockings: 65 % PAD nylon, 35 % elastane LYCRA®

This medical aid does not contain any active substance, including human blood or plasma derivatives, tissues or cells of human origin or derivatives thereof, or tissues or cells of animal origin or derivatives thereof.

Disposal:

Compression medical stockings can be disposed of with normal household waste or in a textile container. Although the product is intended for use on healthy, intact skin, it may become contaminated with body fluids during use. In such a case, dispose of the compression medical aid as infectious material.

The material used is non-toxic unless ignited. Never incinerate used compression medical aids.

Any serious adverse event occurring in connection with the product in question should be reported to the manufacturer and to the competent authority of the Member State where the user and/or patient is located:

Czech Republic	Státní ústav pro kontrolu léčiv, Šrobárova 48, 100 41 Prague 10
Slovak Republic	Štátny ústav pre kontrolu liečiv, Kvetná 1024/11, 821 08 Bratislava
Poland	Urząd Rejestracji Produktów Leczniczych, Wyrobów Medycznych i Produktów Biobójczych, Al. Jerozolimskie 181C, 02-222 Warsaw
Germany	Das Bundesinstitut für Arzneimittel und Medizinprodukte, Kurt Georg Kiesinger Allee 3, 53175 Bonn
Russia	Roszdraznadzor, Address: 4, bld. 1, Slavyanskaya Square, Moscow, 109074
Ukraine	Ministry of Health Ukraine, 7 Hrushevskoho Street, Kiev, 01601
Australia	Therapeutic Goods Administration, 136 Narrabundah Lane, Symonston ACT 2609

List of abbreviations:

KT, CCL – compression class

mmHg – millimeters of mercury column (unit of blood pressure measurement on the leg)

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