Avicenum **ORTHO 360**

(EN) INSTRUCTION MANUAL

Avicenum ORTHO 360 wrist sleeve type 01

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(EN) Dear customer,

In your hands you are now holding a wrist sleeve for supportive post-operative and post-traumatic therapy. You can maximize the benefits of the product by its proper use and maintenance.

Avicenum ORTHO 360 Wrist Sleeve Type 01 is manufactured by the circular knitting technology. The materials used are non-irritating and latex-free. The product achieves the desired therapeutic effect by exerting compression pressure in the wrist area (you must choose the correct size). It is intended for wrist support in post-traumatic and post-operative therapy. The product is most effective during physical strain. The product is also suitable for sports activities. You should remove the product during longer breaks. **Avicenum ORTHO 360 Wrist Sleeve Type 01** is a medical device intended for use both at home and in medical facilities.

Product name:

Avicenum ORTHO 360 Wrist Sleeve Type 01

Product variants:

- double-sided, sizes S – XXL

User profile:

Avicenum ORTHO 360 Wrist Sleeve Type 01 is a medical device intended for use both at home and in medical facilities.

Intended use:

A medical device designed to support the wrist during post-traumatic and post-operative therapy, in case of flabby muscles, inflammation, or injury to wrist ligaments and tendons and to prevent swelling. This product is also suitable for sports activities as a therapeutic or preventive aid.

Indications:

- wrist support in post-operative or post-traumatic therapy
- instability and weakness of wrist ligaments
- sprain or injury of wrist ligaments and tendons
- inflammation and swelling in the wrist area
- carpal tunnel syndrome
- prevention of injury and wrist strengthening (e.g., during sporting activities)

If this is the first time you use a wrist sleeve, consult your doctor or have your medical condition regularly examined.

Contraindications:

- damaged skin
- acute weeping skin manifestations
- lymphatic swelling
- hypersensitivity to materials used
- limb sensitivity disorders (e.g., peripheral neuropathy)

Your doctor should consider the use of the product in view of the known contraindications and possible risks compared to the product's therapeutic benefits.

Target group:

Patients with wrist instability, patients undergoing post-operative or post-traumatic therapy of the wrist, and patients with the above indications. The product is intended for use on healthy, intact skin. The product is intended for all age groups and both sexes. The product may be used in pregnancy without limitation.

Side effects:

Side effects may occur especially when using a wrong size of the product or if the product is used outside strenuous activities. If the product is too tight, there may be pressure sores, impaired blood circulation due to compression of blood vessels, or deterioration of sensitivity due to excessive compression of the nerves. Sensitive and hypersensitive individuals may experience skin irritation (redness, itching, or blistering) at the point of contact of the product with the skin.

How to choose the correct size:

In order to achieve the desired effects, you must select the correct size of the product. Determine the size by measuring the circumference at the points indicated in the figure.

Note:

Any swelling or instability of the wrist should be examined by your doctor. If this is the first time you experience swelling or instability of the wrist, consult your doctor before using the product. Do not use the product if you do not know the cause of the swelling or instability of the wrist. Always consult your doctor before combining the product with other medical devices (e.g., compression arm sleeves).

- use a tape measure
- take the measurements standing up, preferably in the morning right after waking up
- get assistance from another person
- measured values may change during the therapy – always verify the correct size before using a new product

Size	Wrist circumference in cm (point c)
S	13–15
M	15–17
L	17–19
XL	19–21
XXL	21–23

If you are on the border between two sizes, choose the bigger one. If you have difficulties choosing the size of the product, visit a specialized medical supplies store or directly the manufacturer.

If you have difficulties putting on the product, visit a specialized medical supplies store and request a demonstration of the correct method of applying it.

With a correctly selected size:

- you should feel noticeable but not unpleasant pressure at the wrist area
- the product does not slide or roll down
- the product does not choke at any point

Size and manufacturer information can be found on the tag sewn at the top edge of the product.

Method of use:

Put on the sleeve by gradually pulling it over the wrist and poking the thumb through the side hole. Never put on the sleeve by merely pulling the top edge – always put it on gradually. When putting on, the product may not come into contact with sharp objects (e.g., jewellery) or long nails.

Stretch the product evenly without creases for effective distribution of pressure.

If you use body cosmetics, wait about 15 minutes after application before you put on the product.

In the case of visible mechanical damage to the product, stop wearing it and preferably replace it with a new product.

Care and maintenance:

Avicenum ORTHO 360 Wrist Sleeve Type 01 is a medical device intended for repeated use by a single person. The product will maintain its properties and safety for 12 months from the first use under these conditions:

- the product can be washed in a washing machine (maximum washing temperature 30 °C, gentle cycle) in a soapy solution or using a special preparation for washing compressing stockings and sleeves using a wash bag
- do not use fabric softeners
- do not use centrifugation; squeeze the wet product between two towels
- dry the product in a horizontal position away from a direct heat source (radiator, sun)
- do not iron the product
- the maximum number of washing cycles is 120
- store in a dry and dark place, preferably in the original packaging
- do not interfere mechanically in the product, do not repair places damaged by use, do not cut the edges of the product
- the product may not come into contact with organic solvents or bleaching agents

The product is safe to use for 60 months from the date of manufacture. This date is included in the product's batch number in the format YY/MM-1234567 where YY is the year and MM the month of manufacture.

Care symbols:

Do not use fabric softeners!

Prohibited use:

Stop using the product immediately if:

- If it shows signs of mechanical damage (replace the product)
- Your health has deteriorated which could be associated with the use of the product – consult your doctor at once

Composition:

60 % PAD Nylon, 40 % elasthan LYCRA®

This medical device does not contain an active substance, including human blood or plasma derivatives, tissues or cells of human origin or derivatives thereof, tissues or cells of animal origin or derivatives thereof.

Disposal:

The product can be disposed of with normal household waste or in a textile waste container. Although the product is intended for use on healthy and intact skin, it may be contaminated with body fluids during use. In that case, dispose of the product as an infectious material.

The material used is not toxic unless burned. Do not burn the used sleeve.

Any serious adverse event that occurred in connection with the use of this product should be reported to the manufacturer and the competent public authorities of the user's country of residence.

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Україна	Ministry of Health Ukraine, 7 Hrushevskoho Street, Kyiv, 01601
Australia	The Therapeutic Goods Administration, 136 Narrabundah Lane, Symonston ACT 2609

Manufacturer:

ARIES, a.s., Studenec 309, 512 33 Studenec, Czech Republic
Last reviewed: 21.06.2024

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