# *Avicenum*

## **COTTON 360**

### **PHLEBO**

#### **(GB) Instructions for Use**

Avicenum PHLEBO 360 COTTON Calf-Length Stockings  
Avicenum PHLEBO 360 COTTON Thigh-High Stockings

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Avicenum PHLEBO 360 COTTON Wadenstrümpfe  
Avicenum PHLEBO 360 COTTON Schenkelstrümpfe

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**(GB)**

## **INSTRUCTION MANUAL**

Dear customer, In your hands you are now holding compressive medical stockings that will help you in the treatment of venous or lymphatic diseases. You can maximize the benefits of the product by its proper use and maintenance.

**Avicenum PHLEBO 360 COTTON** are compressive medical stockings of compression class II (KT II, CCL 2). The product was manufactured by the circular knitting technology. The materials used for the manufacture of the product are non-irritating and latex-free. Thanks to the high content of natural cotton (60%), the product is particularly suitable for persons with sensitive skin or allergy. The mechanism behind the effects of the product is based on the action of precisely defined graduated pressure on the leg, which reduces from the ankle up to the heart. The pressure applied in the ankle area (see the size table) is the greatest and reaches 23–32 mmHg. In order to achieve the desired effects, you must select the correct size of the product. The product is intended for physical exercise. You should take off the product for extended periods of rest (unless your doctor determines otherwise). **Avicenum PHLEBO 360 COTTON** is a medical device intended for use both at home and in medical facilities.

**Product name:**

### **Avicenum PHLEBO 360 COTTON Calf-Length Stockings**

**Product variants:**

- open toe, KT II, sizes S normal – XXL long
- closed toe, KT II, sizes S normal – XXL long

### **Avicenum PHLEBO 360 COTTON Thigh-High Stockings**

**Product variants:**

- lace, open toe, KT II, sizes S normal – XXL long
- lace, closed toe, II KT II, sizes S normal – XXL long
- hem, open toe, KT II, sizes S normal – XXL long
- hem, closed toe, KT II, sizes S normal – XXL long

**Product accessories:**

Each package of **Avicenum PHLEBO 360 COTTON** in all variants includes a textile aid for putting on Avicenum compression hosiery. The aid is not a medical device and only serves to make it easier to put on the product. The correct use of the aid is described below, together with the instructions for the proper wearing of the compression hosiery.

**User profile:**

**Avicenum PHLEBO 360 COTTON** is a medical device intended for use both at home and in medical facilities.

**Intended use:**

The product is intended for the compression therapy of venous conditions and swelling of the legs of various origin.

**Indications:**

- chronic venous insufficiency – subjective and objective manifestations, prevention
- primary and secondary varices of the legs
- deep venous thrombosis
- superficial venous thrombophlebitis
- swelling associated with post-thrombotic syndrome after deep venous phlebothrombosis
- swelling and varices during pregnancy
- conditions after phlebothrombosis and thrombophlebitis of the legs
- prevention of deep venous phlebothrombosis in surgery and associated thromboembolic diseases (especially surgery, orthopaedics, gynaecology)
- swelling of the legs (post-operative and post-traumatic, lymphatic, overweight swelling, idiopathic swelling of the legs)
- compression support after sclerotherapy of leg varices
- compression support after surgical extirpation of varices (standard surgery, radiofrequency ablation, or laser ablation)
- venous compression syndromes (Cockett's syndrome – compression of the left iliac vein by the right iliac artery, compression of the popliteal vein, surrounding compression of veins, e.g., by tumour, cyst, or surrounding haematoma)
- venous malformations (Klippel-Trenaunay syndrome, Parkes Weber syndrome, Maffucci syndrome)
- lymphedema (primary and secondary)
- lipoedema
- leg ulcer therapy
- travel thrombosis (economy class syndrome)

If this is the first time you use compression medical hosiery, consult your doctor to select the proper compression class of the product or have your medical condition regularly examined.

**Contraindications:**

- Acute weeping skin manifestations
- acute limb ischemia
- critical limb ischemia
- decompensated cardiac insufficiency
- phlegmasia coerulea dolens, phlegmasia alba dolens
- Hypersensitivity to materials used
- limb sensitivity disorders (e.g., peripheral neuropathy in diabetes mellitus)

Your doctor should consider the use of the product in view of the known contraindications and possible risks compared to the product's therapeutic benefits.

**Target group:**

Patients with primary varices, patients after venous thrombosis, and patients with the above indications. The product may be used only on healthy and intact skin. The product is intended for all age groups and both sexes. The product may be used in pregnancy without limitation.

**Side effects:**

Side effects may occur especially when using a wrong size of the product or if the product is used outside strenuous activities. This may cause pressure sores, impaired blood circulation due to compression of blood vessels, or deterioration of sensitivity due to excessive compression of the nerves. Despite the high content of cotton, sensitive and hypersensitive individuals may experience skin irritation (redness, itching, or blistering) at the place of contact of the product materials with the skin.

**How to choose the correct size:**

In order to achieve the desired effects, you must select the correct size of the product. Determine the size according to the size table by measuring the circumference at the points indicated in the figure.

**Note:**

Any leg swelling should be examined by your doctor. If this is the first time you experience the swelling, consult your doctor before using the product. Do not use the product if you have not determined the cause of the swelling. Always consult your doctor before combining the product with other medical devices (e.g., sleeves or braces).

- Use a tape measure
- Take the measurements standing up, preferably in the morning right after waking up – the legs may swell during the day
- Get assistance from another person
- Measured values may change during the therapy – always verify the correct size before using a new product


**SIZE TABLE:**

SIZE	a	b	b1	c	d	e	f	g
S	19 – 25	20 – 23	25 – 30	32 – 38	31 – 37	33 – 39	41 – 49	47 – 57
S+	19 – 25	20 – 23	25 – 30	32 – 38	31 – 37	33 – 39	49 – 57	57 – 67
M	22 – 28	23 – 26	28 – 33	35 – 41	34 – 40	36 – 42	45 – 53	52 – 62
M+	22 – 28	23 – 26	28 – 33	35 – 41	34 – 40	36 – 42	53 – 61	62 – 72
L	25 – 31	26 – 29	31 – 36	38 – 44	37 – 43	39 – 45	49 – 57	57 – 67
L+	25 – 31	26 – 29	31 – 36	38 – 44	37 – 43	39 – 45	57 – 65	67 – 77
XL	28 – 34	29 – 32	34 – 39	41 – 47	40 – 46	42 – 48	53 – 61	62 – 72
XL+	28 – 34	29 – 32	34 – 39	41 – 47	40 – 46	42 – 48	61 – 69	72 – 82
XXL	31 – 37	32 – 35	37 – 42	44 – 50	43 – 49	45 – 51	57 – 65	67 – 77

LENGTH		
	normal	long
AD	34 – 38	38 – 43
AG	62 – 71	71 – 82

- a foot circumference
- b circumference above the ankle
- b1 circumference 10 cm above the ankle
- c calf circumference
- d circumference below the knee
- e circumference above the knee
- f thigh circumference (in the middle of the thigh)
- g thigh circumference (5 cm under the crotch)
- h hip circumference
- t waist circumference

- AD calf-length stockings
- AG thigh-high stockings

If the measured values do not correspond to any of the sizes in the table (e.g., point b corresponds to size S but point d is already size L), visit a specialized medical supplies store or the manufacturer. It could be that you need a custom-made product.

If the measured values differ significantly for the right and the left leg, visit a specialized medical supplies store or directly the manufacturer. It could be that you need a custom-made product.

If you are on the border between two sizes, choose the bigger one. If you have difficulties choosing the size of the product, visit a specialized medical supplies store or directly the manufacturer.

If you have difficulties putting on the product, visit a specialized medical supplies store and request a demonstration of the correct method of applying it. You can also use one of the aids for putting on the product.

If the hem or lace does not hold on the leg (in thigh-high stockings), clean the silicone straps or studs with an alcohol solution or shave the hair at the point of contact of the skin with the silicone straps or studs.

**With a correctly selected size:**

- you should feel noticeable but not unpleasant pressure over the entire length of the product
- the product does not slide nor roll down below the knee or on the instep when walking
- the lace or hem does not choke (in thigh-high stockings)

Information about the compression class, size, and manufacturer can be found on the tag sewn at the top of each product. If the tag is removed over the course of using the product, its size is still indicated on the sole of the product.

#### **Method of use:**

Put on the product by gradually pulling it over the leg (do not pull only the top edge). It is preferable to put on the product in the morning before any physical strain or formation of swelling. To ease up the process, use the supplied textile aid. The method of use of the aid is illustrated on the card supplied with the product. When putting on, the product may not come into contact with sharp objects (e.g., jewellery) or long nails. Intensive care for your feet (no rough skin on the heels) will decrease the risk of damaging the product and increase its life. Stretch the product evenly without creases for effective distribution of pressure. If you use body cosmetics, wait about 15 minutes after application before you put on the product. In the case of visible mechanical damage to the product, stop wearing it and preferably replace it with a new product.

#### **Care and maintenance:**

**Avicenum PHLEBO 360 COTTON** is intended for repeated use by a single person. The product will maintain its properties and safety for six months from the first use under these conditions:

- The product can be washed in a washing machine (maximum washing temperature 30 °C, gentle cycle) in a soapy solution or using a special preparation for washing compression hosiery using a protective container (wash bag)
- Do not use fabric softeners
- Do not use centrifugation; squeeze the wet product between two towels
- Dry the product in a horizontal position away from a direct heat source (radiator, sun)
- Do not iron the product
- The maximum number of washing cycles is 120
- Store in a dry and dark place, preferably in the original packaging
- Do not interfere mechanically in the product, do not repair places damaged by use, and do not cut the hem or lace of the product
- You can restore the adhesive properties of the silicone strips on the lace or of the silicone studs on the hem by applying an alcohol solution (in thigh-high stockings)
- The product may not come into contact with organic solvents or bleaching agents

The product is safe to use for 60 months from the date of manufacture. This date is included in the product's batch number in the format YY/MM-1234567 where YY is the year and MM the month of manufacture.

#### **Care symbols:**

*Do not use fabric softeners*

**Prohibited use:**

Stop using the product immediately if:

- if it shows signs of mechanical damage (replace the product)
- your health has deteriorated which could be associated with the use of the product – consult your doctor

**Composition:**

Calf-length stockings: 60% cotton, 5% PAD Nylon, 35% elasthan LYCRA®

Thigh-high stockings: 60% cotton, 5% PAD Nylon, 35% elasthan LYCRA®, silicone

This medical device does not contain an active substance, including human blood or plasma derivatives, tissues or cells of human origin or derivatives thereof, tissues or cells of animal origin or derivatives thereof.

**Disposal:**

The product can be disposed of with normal household waste or in a textile waste container. Although the product is intended for use on healthy and intact skin, it may be contaminated with body fluids during use. In that case, dispose of the product as an infectious material.

**The material used is not toxic unless burned. Never dispose of the product by burning it. Any serious adverse event that occurred in connection with the use of this product should be reported to the manufacturer and the competent public authorities of the user's country of residence.**

Czech Republic	Státní ústav pro kontrolu léčiv, Šrobárova 48, 100 41 Prague 10
Slovak Republic	Štátny ústav pre kontrolu liečiv, Kvetná 1024/11, 821 08 Bratislava
Poland	Urząd Rejestracji Produktów Leczniczych, Wyrobów Medycznych i Produktów Biobójczych, Al. Jerozolimskie 181C, 02-222 Warsaw
Germany	Das Bundesinstitut für Arzneimittel und Medizinprodukte, Kurt Georg Kiesinger Allee 3, 53175 Bonn
Russia	Roszdraznadzor, Address: 4, bld. 1, Slavyanskaya Square, Moscow, 109074
Ukraine	Ministry of Health Ukraine, 7 Hrushevskoho Street, Kiev, 01601
Australia	Therapeutic Goods Administration, 136 Narrabundah Lane, Symonston ACT 2609

**Abbreviations:**

KT, CCL – compression class

mmHg – millimetres of a column of mercury (unit for measuring blood pressure in the legs)

**Manufacturer:**

ARIES, a.s., Studenec 309, 512 33 Studenec  
Last reviewed: 01.03.2022

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**License Partner  
Certificate**

We license only to:  
**Aries A.S.**

For continued all necessary steps that the respective activities have no harmful impact on the environment for the category "Textiles". Further GOTS and Fair for a certified license requires certain company.

The company has been approved as a license partner and has the right to designate its respective products with the GOTS and Fair for a certified license.

INTERNATIONAL COOPERATION OF TEXTILE RESEARCH ASSOCIATION

**OEKO-TEX®**  
CERTIFICATION IN TEXTILES

## CERTIFICATE

The company

**ARIES, a.s.**  
51233 Studenec u Brna, CZECH REPUBLIC

is certified according to standards EN ISO 9001:2015 and ISO 14001:2015 and to GOTS 100, Fair for a certified license, according to the year 1999-2019.

**OEKO-TEX®**  
STANDARD 100  
by Hohenheim University  
of Applied Sciences  
Textile Technology  
and Fashion Design  
Institute

for the following activities:

Manufacture and supplying goods (textiles, non-textiles, conditioned sport products, quilt and leather items) and trade in the EU, UK, US, and the Middle East. Also, some services and trade of imported goods supplied by OEKO-TEX® producer that includes goods, printed clothing in total and all OEKO-TEX® Standard 100. Products of the respective work activities (textiles) and a GOTS 100, Fair for a certified license product that is also not used and materials and accessories (textiles) of OEKO-TEX® Standard 100, Fair for a certified license. The goods are not used in the EU.

For certain other activities (products of the EU of OEKO-TEX® for use in non-textiles, other items, etc.) and for the goods that are not used in the EU.

In case of any doubt, see the label's certificate information according to GOTS 100. License is valid for the period 01.01.2019 to 31.12.2019 and is a certificate with process that includes all changes and updates. The certificate is valid for 10 years.

**The certificate 00446 is valid until 31.12.2019**

Date: 14.06.2017

OEKO-TEX® approved for all long-term business

**Aries**  
MEDICAL PRODUCTS

ARIES, a.s., 512 33 Studenec 309  
MADE IN CZECH REPUBLIC  
ISO 9001:2015

[www.avicenum.eu](http://www.avicenum.eu)

SD 11-023 (R03)